

FIFTY50[®]

Half Our Profits Fund Diabetes Research

Glucose Meter 2.0

Next Generation Blood Glucose Monitoring

No-Coding Convenience and Reliability

In a feature-rich,
easy-to-use meter.



Actual Size

Features

- ▶ No-Coding required
- ▶ Simple, one step operation
- ▶ Fast, 5 second test time
- ▶ Small, 0.5 μ l sample size
- ▶ Large, easy-to-read display
- ▶ Alternate site testing
- ▶ Pre-/Post-meal flagging
- ▶ Pre-/Post-meal averaging
- ▶ Three user configurable alarms.
- ▶ USB data port

Product Specifications

Reported result range	20-600mg/dL (1.11 ~ 33.3mmol/L)
Sample size	Minimum 0.5 μ l
Test time	5 seconds
Sample type	Fresh capillary whole blood
Calibration	Plasma equivalent
Assay method	Electrochemical
Battery life	1,000 tests
Power	Two 3.0 V lithium batteries (type CR2032)
Memory	250 test results
Size	3.66 x 1.85 x 0.59 in. (93 x 47 x 15 mm)
Weight	50g (with battery)

Operating Ranges

Temperature	50-104° F (10-40° C)
Relative humidity	10 - 90%
Hematocrit	20-60%

The FIFTY 50 Meter 2.0 makes testing quick and easy with results you can trust.

Fast and Gentle Testing

- ▶ Small 0.5 μ l blood sample for virtually pain-free testing
- ▶ Five second test time – no additional time needed for coding
- ▶ Approved for alternate site testing



1 Insert test strip



2 Apply blood sample



3 Result in 5 seconds


Accurate and Reliable Results

- ▶ Proven accuracy
- ▶ Reliable automatic code-reading mechanism
- ▶ No interference with Maltose, Lactose, Galactose, or Xylose

Feature rich, yet simple to use

- ▶ Pre-/Post-meal flagging
- ▶ Pre-/Post-meal averaging
- ▶ Multiple programmable reminders can be set with the press of a button

Clinical Evaluation Results



Compact, Ergonomic Design

